## **✓** HENRY SCHEIN®

# Micro-motor L50



Copyright © 2009 by Henry Schein Inc.
All rights reserved. No part of these instructions may be reproduced or copied in any form by any means – graphic, electronic or mechanical, including photocopying, typing or information retrieval systems – without written permission of Henry Schein Inc. Todos los derechos reservados. Ninguna parte de estas instrucciones puede ser reproducida o copiada de ninguna manera y por ningún medio - gráfico, electrónico o mecánico, incluyendo el fotocopiado, tipeo o los sistemas de recuperación de información – sin la autorización expresa por escrito de Henry Schein Inc.

Tous les droits réservés. Aucune partie de ces instructions ne peut être reproduite ou copiée par aucun moyen – graphique, électronique ou mécanique, y compris la photocopie, la copie à main ou les systèmes de récupération de l'information – sans l'autorisation par écrit de Henry Schein Inc.

Alle Rechte vorbehalten. Dieses Benutzerhandbuch darf ohne schriftliche Genehmigung von Henry Schein Inc weder ganz noch in Auszügen in jedweder Form- grafisch, elektronisch oder mechanisch- reproduziert oder kopiert werden, was Fotokopien, Abschriften oder das Weiterleiten über andere Kommunikationsmittel einschließt.

Tutti i diritti riservati. Nessuna parte di queste istruzioni può essere riprodotta o copiata in nessuna forma e con nessun mezzo - grafico, elettronico o meccanico, includendo la fotocopiatura, la battitura o sistemi di recupero delle informazioni - senza il permesso scritto della Henry Schein Inc.

Alle rechten voorbehouden. Deze instructies of een deel ervan mogen niet worden gereproduceerd of gekopieerd in welke vorm dan ook d.m.v. grafische, elektronische of mechanische (inclusief fotokopiëren, typen of informatieherstel) systemen zonder de schriftelijke toestemming van Henry Schein Inc.

✓ HENRY SCHEIN®

900-2963

# MICRO-MOTOR L50

•Brushless micro motor • Knee or foot control • 50.000rpm

# **OPERATING MANUAL**

# MICROMOTOR SET L50

• Micromotor escobillas • Control con rodilla o pié • 50.000 rpm

## **INSTRUCCIONES DE USO**

## MICRO-MOTEUR SET L50

- Micro-moteur sans charbon
- Contrôle au genou ou pied 50.000 tr/mn

## **MODE D'EMPLOI**

# HS-Mikromotor Set L50

- Kollektorloser Mikromotor Knie- oder Fußschalter
- 50.000 UpM

## **BENUTZERHANDBUCH**



## MICROMOTORE L50

- Micromotore da laboratorio
- Controllo a pedale o del ginocchio 50.000 rpm

## MANUALE D'USO

**MADE IN KOREA** Fabricado en Corea Fabriqué en Corée Hergestellt in Korea Prodotto in Corea





Distribuido por: Distribué par:

Vertrieb durch: Distribuito da

Micro

HENRY SCHEIN INC. Melville, NY 11747 USA

EC REP

Henry Schein U.K. Holdings Ltd. Gillingham ME8 0SB U.K.

www.henryscheinbrand.com



## ✓ HENRY SCHEIN®

# Micro-Motor L50

Thank you for purchasing our product.

Proper maintenance and repair extends the life of the product.

Keep the instructions for use in a safe place for your reference.

Accidents caused by not following the instructions for use will be deemed not our responsibility.

Product appearance or specification may be changed without notice. The users are responsible for the operation and maintenance of the product.

In the event of a problem please state model name and full details to enable us to provide prompt attention.

## WARRANTY

Handpiece - 1 year

Control box & Foot Pedal - 1 year

Warranty does not include misuse or normal wear of bearings, carbon brush, collet chuck, spindle and armature.

### Other cases

- Act of God (fire, sea wind or water, flood etc.).
- Replacing by the normal wear of the consumable parts (ball bearings, collet chuck, spindle, or armature etc.).

## SERVICE CHARGE

In cases of fault or malfunction due to misuse by the user, for example:

- Using with wrong Input Voltage (AC Voltage).
- Misuse or improper use.
- Dropping during setting or moving.
- Using parts or accessories not recommended by manufacturer.
- Malfunction if repaired by other than manufacturer recommended repairers.

## ✓ HENRY SCHEIN®

# PRODUCT SPECIFICATION

## **MICRO-MOTOR L50 (KNEE TYPE)**

INPUT POWER SOURCE AC220-230V 50/60HZ MINIMUN LOAD CURRENT 2A

WEIGHT 3Kg

DIMENSIONS (mm) W 95 - L 270 - H 230 100 5A - 220/230 2A

FUSE


### **HANDPIECE**

NON LOAD RPM 50,000 RPM OUTPUT 230W 7.8 N.cm MAX. TORQUE

WEIGHT 254a

DIMENSIONS (mm) L 171 - L1 78 - L2 93

## PART NAME



- CONTROL BOX
- KNEE CONTROLLER
- **VOLUME DIAL**
- FORWARD-REVERSE BUTTON
- **AUTO-CRUISE BUTTON**
- HAND(KNEE)-FOOT BUTTON
- RPM DISPLAY BOX
- POWER SWITCH
- 9. MOTOR CONNECTION
- 10. FOOT PEDAL CONNECTION

- 11. POWER CONNECTION
- 12. FUSE BOX
- 13. DEVICE FOR ATTACHMENT (Wall mounting)
- 14. HANDPIECE
- 15. POWER CORD
- 16. HANDPIECE STAND
- 17. WRENCH
- 18. FOOT SPEED CONTROL PEDAL
- 19. BOLTS FOR DEVICE
- \*Please check all items are included before use.

<sup>\*</sup>Design and specification of the product may be subject to be change.

# SETTING & OPERATION


Please attach the knee controller as per the image below. Avoid any excess protrusion of the bolts.

Ensure that all the bolts are tightened.

## POWER ON (I)

Turn on the power switch.


The unit will beep, the switch will illuminate green and the readout will show on the RPM display box.



### POWER OFF (O)

Unit will shut down.

- 1. Check the power is **OFF**.
  - A Connect the power cord.
  - **B** Connect the foot control pedal to the foot pedal connector.
  - **C** Connect the handpiece to the handpiece connector.
- 2. Set the speed control dial on the control box to '0'.
- 3. Check the collet chuck is locked by checking the fitting or try to pull out the bur.
- 4. Turn on the power switch.
- **5.** Select Hand(Knee)/Foot control by pressing the button . Each press will change the selection.
- **6.** Select direction by pressing the forward/reverse button. Each press will change the selection.



**RPM Chart** 

### In case of overload

Rpm display box shows 00 and control box stops after 5 seconds displaying Error 1 with beeping sound. To release, push motor button when working by hand or lift off foot pedal.

# MAINTENANCE

- 1. Clean the handpiece surroundings.
- **2.** Wipe handpiece with disinfection wipes of isopropyl alcohol, avoid moisture inside the handpiece. Do not apply the oil to the unit. This contaminates the grease inside of the bearing. Do not clean with water.
- **3.** Air cleaning is acceptable, but care should be taken as too much air pressure may affect to the handpiece.
- **4. Bur insertion and removal:** Turn the chuck until its clicks either in the open direction and the bur can be removed or in the locking direction and the bur is fixed.

Do not overturn or this will result in the handpiece becoming disassembled.

## TROUBLESHOOTING

	PROBLEM	CHECK	CAUSE	SOLUTION
	The motor does not run	The power lamp remains off, when the power switch is on.	Control box may have a problem.	Check the fuse. Check the power switch.
		The handpiece does not run, when the hand/foot switch is on and power lamp is on.	Check the power cord connection.	Connect the power cord properly.
		The handpiece does not run, when the power lamp is on.	Check if the chuck is open. (E1 error)	Lock the chuck.
			Check the connection between the CB jack and HP cord.	Connect the Handpiece properly.
			Check if the cord is broken down or disconnected.	Consult the supplier.
		The RPM does not increase.	The control box may have a problem.	Consult the supplier.

If the product does not work after trying the above, please contact to the supplier or manufacturer.

4

MICRO-MOTOR L50

MICRO-MOTOR L50

# MAINTENANCE

1. Memory:

Saves the previous selection of direction or hand/foot operation. Exercise care on next use of the unit.

2. Do not use the control box in humid or damp environment.

3. Transport and storage conditions:

Temperature: 10 ~ 40°C Humidity: 30 ~ 75% RH

Atmospheric Pressure: 700 hPa to 1060 hPa

- 4. Use the handpiece and pedal only with the control box provided.
- 5. Check for any problems, including vibration of bur, at low speed before selecting high speed.
- 6. Operate the speed control dial control box gently.
- 7. Turn off the power or set the RPM is to zero and then change the bur.
- 8. Do not turn on the power while the collet chuck of handpiece open. (E1 error).

### 9. Maximum operating temperature of component parts:

Head of the Handpiece - maximum temperature 39°C at ambient 25°C or 54°C at ambient 40°C Therefore use carefully!!

Recommended RPM according to the bur diameter

Please keep the RPM according to below diagram to avoid accidents.

(The allowable RPM could be different for the bur material and type.)

### **10.** Bur

- Bur dimension: 2.35mm (ISO7785-2 TYPE2)
- · Minimum fitting length of bur: 10mm
- · After insertion into the handpiece and locked in position, the maximum projection length of bur should be at least 35mm.

Use ISO standard burs and avoid using rusty and bent burs which will cause excessive wear

Operate at speed recommended by bur manufacturer.

- 11. Do not to drop the handpiece as this may result in damage or breakage. If dropped, check the bur is not bent or that there is no excessive heat is generated while operating.
- 12. If any problems occur in handpiece, have handpiece checked by a designated repairer before use.
- 13. ENSURE POWER SWITCH OF THE UNIT IS OFF BEFORE CONNECTING TO THE MAINS.
- 14. A circuit breaker is provided to protect the motor. It activates when there is an overload.
- 15. Do not open the chuck while operating.
- **16.** Keep the handpiece with a bur within the collet chuck after operation.
- 17. Do not put the handpiece down before it has completely stopped.
- **18.** Remove the power cord from the supply before removing the handpiece.
- 19. Use the specified voltage.
- 20. Only use accessories suitable for this product.
- 21. For recommended speeds in various applications, follow the manufacturer's instructions.

### **CHANGING the FUSE**

The fuse is placed in the small box next to the power connector at the back of the unit. Remove the fuse box and change the fuse.

# Error Code

In case of any overload, breakage, and motor coil burning, the motor stops and control box displays an error code on the rpm display box. Diagnosis and corrective action can be checked from this:

CODE	CAUSE	RELEASING	ACTION
E1	Setting error: Caused by abnormal of motor initial operation.	Push the motor button, knee or foot controller.	Check the power cord and motor connector are correctly connected. Check the chuck is not open. Turn the power on again. If operating correctly the error code will not display. If error persists contact the supplier/manufacturer.
E2	Over current error: Occurs in case of excessive use of handpiece on overload. Incorrect assembly. Bearings and motor wear.		Switch off the control box and restart after 10 minutes. If error persists contact the supplier/manufacturer.
E3	Control box error.		Contact the supplier/manufacturer.

### **CHANGE COLLET CHUCK**

✓ HENRY SCHEIN®

- 1. Set and fix the wrench turning the bur.
- 2. Turn in open direction
- 3. 4. Turn the chuck to the counter-clock wise by hand or if difficult use wrench.
- 5. Change the chuck.

\* Fix the bur and chuck turning to clock wise. And then turn again to counter clock wise a bit and set the tension to remove the bur. Clean and assemble the chuck very carefully.

If your are not a skilled engineer, do not disassemble the handle and motor, It may cause a breakage,

If disassembled, the joint of handle part and the shaft of the motor part have to be connected properly.

If it does not fit, it may cause vibration, noise, heat and bearing damage.

\* If the chuck is not open, please turn it with the wrench.

For your safety, please note the following:

Do not use with damaged or faulty power cord, plug, and outlet, this may cause the danger of fire/ electric shock.

When unplugging the power cord, hold the plug and do not touch with wet hand - this may result in an electric shock

Do not bend the power cord excessively or crush with heavy objects as this may cause fire/electric shock.

Do not place near heat, candles or cigarettes

Close supervision is needed if this product is used by or near the very young, the elderly or invalids.

Do not leave children unattended in any room where this product has been installed.

- · Type of protection against electric shock: Class I equipment
- Degree of protection against electric shock: Type B applied parts
- · Classification according to the degree of protection against ingress of water as detailed in the current edition of IEC 60529: Equipment: IPX0, ordinary equipment

Foot Switch: IPX1

- · This equipment is not suitable for use in the presence of flammable anesthetics or oxygen.
- · Mode of operation: continuous operation

This equipment has been tested and found to comply with the limits for medical devices in IEC 60601-1-2:1994. These limits are designed to provide reasonable protection against harmful interference in a typical medical installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to other devices in the vicinity. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to other devices, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- · Reorient or relocate the receiving device.
- · Increase the separation between the equipment.
- Connect the equipment into an outlet on a different circuit to where the other device(s) are connected.
- · Consult the manufacturer or field service technician for help.

Use the grounding conductor/clip. Breakage or electric leak may cause the danger of electric shock. For grounding, inquiry to electrician or to main office.

Do not use the product in a humid or damp location - this may cause danger of electric shock and fire.

Do not drop or damage the handpiece. In case of malfunction, dropping or immersing in water (this may cause vibration of the bur, overheating, or damage of the bearings.) Do not use - contact the manufacturer,

- \* Only use the product as specified in the instructions of the manufacturer.
- \* Injury or property damage may result if these instructions are not correctly followed.

The manufacturer in not responsible for any accidents resulting from the instructions not being followed.