



blast

Let's get flossing

A convenient and effective way to get your patients to experience the benefits of regular flossing in their daily routines



Available from



Visit henryschein.co.nz for more details

*vs string floss and vs brushing alone.



blast water flosser

The fun and easy way for your patients to floss




Available from



Only **1 in 4**
 Australians **floss daily***

Meet **blast**

Help your patients create a
 better flossing routine



Blasts away plaque
 between teeth and
 along the gums



Gentle cleaning
 and helps improve
 gum health



Easy for patients to
 clean around braces
 and dental work

What makes the Colgate Blast Water Flosser different?

Its portable and uniquely patented design, combined with easy to use controls will encourage your patients to clean their in-betweens and achieve a cleaner feeling than brushing alone.

The Blast Water Flosser has been designed to be simple and effective for people who are new to water flossing with its easy power flossing that's good and gentle for gums.

*ADA 2020 Consumer Survey - a national survey of 25,000 adults over 18 years of age which found that three quarters of Australian adults rarely or never floss or clean between their teeth.


SERIES 1

The fun and easy way to floss


 3 pressure modes
 tailored to you


 Built-in nozzle
 compartment


 Up to 60 secs
 of water capacity^


 Easy to
 fill and clean



 Water resistant


 Rechargeable
 battery



SERIES 2

Your oral health routine upgraded


 Interactive
 LED display


 Large 180mL
 water tank



 Compact
 portable design

Personalise your flossing with **3 specialised nozzles**



Regular

Blast away between teeth
 to combat tooth decay
 for a dentist clean feel



Interdental

Creates a unique
 stream pulse water flow
 to blast away plaque



Orthodontic

An easy clean for those
 hard-to-reach areas
 around braces


Choose from **5 cleaning modes**



Gentle



Regular



Deep Clean



Massage



Spot Spray



^When tank is full and used on regular mode.