



COLGATE® PULSE

SERIES 1

Tracks how long and often you brush to improve your brushing routine

2 modes, 2 minute timer and 2 week battery life†



CG-61030574

CG-61030573

SERIES 2

Tracks how long, often & well you brush, coaching you to a healthier smile

3 modes, pressure sensor, 2 minute timer, 2 week battery life† and travel case



CG-61030577

CG-61030576



REPLACEMENT BRUSH HEADS

Not available at Henry Schein. Only available across select national retailers.



SENSITIVE

Tough on plaque, gentle on gums



WHITENING

Whiter teeth in 1 week*



DEEP CLEAN

Removes 5x more plaque†



PULSE

SERIES 2

SMART COACHING FOR A SUPERIOR CLEAN*



Available from



*at gum margin after brushing for 6 weeks with vs without Colgate Connect app

Available from



Visit henryschein.co.nz for more details

†when used twice daily

*vs. an ordinary manual toothbrush


†at the gum margin after 4 weeks vs a manual flat-trim toothbrush

MOST PEOPLE MISS OVER **40%** OF THEIR TEETH WHEN THEY BRUSH*


GET CONNECTED AND BRUSH BETTER



GUIDED BRUSHING
16 zone detection with speed and angle feedback (Series 2 only)



SMILE REWARDS
Get into a healthy brushing routine and earn points#



GET REMINDED
Engaging daily in-app content and reminder notifications

FEATURES



3 BRUSHING MODES

A simple switch lets you select the right mode for you: Gentle Clean, Daily Clean, Deep Clean (Series 2 only)



CELEBRATE BRUSHING STREAKS

The Smile Light changes colours to indicate brushing achievements and to encourage more consistency



PRESSURE SENSOR

Uses smart technology to detect over-brushing, red LED lights alert you when brushing too hard (Series 2 only)



2 MINUTE TIMER

White LED lights help guide your pace for a complete 2 minute brushing routine



RECHARGEABLE BATTERY

2 week battery life when used twice daily

REMOVE **50% MORE PLAQUE**** (SERIES 2 ONLY)

*Based on tooth surface coverage data collected by the Colgate Connect app when brushing with a smart electric toothbrush.
**at gum margin after brushing for 6 weeks with vs without Colgate Connect app.

#See colgate.com/apprewards for more information on how you could get rewarded