



EDGEENDO®

# EDGE V-FILE<sub>HT</sub>

EdgeV-File HT™ features our Heat-Treated™ FireWire NiTi and is identical to V-Taper™ 2H. Alleviate the stress of a root canal procedure with our super strong and flexible EdgeV-FileHT™ system! EdgeV-Files HT™ system includes tip sizes 20, 22, 25, 30, 35, 40, 45 and 50. They have a changing variable taper along the length of the file.



**\$62.10**  
per 6 pack

 HENRY SCHEIN®  
DENTAL |  Rely on Us™

To order phone 0800 808 855, between 8.00am and 5.30pm [www.henryschein.co.nz](http://www.henryschein.co.nz)

ISO Size/Taper	21mm	25mm	31mm
20/.06	EN-VF200621	EN-VF200625	EN-VF200631
25/.06	EN-VF250621	EN-VF250625	EN-VF250631
30/.06	EN-VF300621	EN-VF300625	EN-VF300631
35/.06	EN-VF350621	EN-VF350625	EN-VF350631
40/.06	EN-VF400621	EN-VF400625	EN-VF400631
45/.06	EN-VF450621	EN-VF450625	EN-VF450631
50/.06	EN-VF500621	EN-VF500625	EN-VF500631
22/.07	EN-VF220721	EN-VF220725	EN-VF220731
25/.08	EN-VF250821	EN-VF250825	EN-VF250831
30/.08	EN-VF300821	EN-VF300825	EN-VF300831

### GlidePath Files

ISO Size/Taper	17mm	21mm	25mm	29mm	31mm
13/.03		EN-VFGP130321	EN-VFGP130325		EN-VFGP130331
14/.03		EN-VFGP140321	EN-VFGP140325	EN-VFGP140329	
17/.04		EN-VFGP170421	EN-VFGP170425		EN-VFGP170431
25/.09	EN-VFGP250917				

### Standard Gutta Percha & paper points

ISO Size/Taper	Gutta Percha Point \$18.20 pk 60	Paper Point \$16.70 pk 60
15/.04	EN-GP1504	EN-PP1504
15/.06	EN-GP1506	EN-PP1506
20/.04	EN-GP2004	EN-PP2004
20/.06	EN-GP2006	EN-PP2006
25/.04	EN-GP2504	EN-PP2504
25/.06	EN-GP2506	EN-PP2506
30/.04	EN-GP3004	EN-PP3004
30/.06	EN-GP3006	EN-PP3006
35/.04	EN-GP3504	EN-PP3504
35/.06	EN-GP3506	EN-PP3506
40/.04	EN-GP4004	EN-PP4004
40/.06	EN-GP4006	EN-PP4006
45/.04	EN-GP4504	EN-PP4504
45/.06	EN-GP4506	EN-PP4506
50/.04	EN-GP5004	EN-PP5004
50/.06	EN-GP5006	EN-PP5006



**EDGEENDO**®

**HENRY SCHEIN**®  
DENTAL

**Relay on Us**

To order phone 0800 808 855, between 8.00am and 5.30pm

[www.henryschein.co.nz](http://www.henryschein.co.nz)