
Bendable Periosteal Elevator for Surgical Treatment

KimTrac







Shape Memory

After being bent, NiTi models can be returned to its original shape at 60°C or higher

Serrated Edges

Stably secures the cortical bone plate, regardless of its shape, due to the saw-toothed wheel at the end of the KimTrac.

Various Angles

Consists of three types of angles that make contact with the bony plate. It adheres well to both anterior and posterior locations by adjusting the left-right angle.

KimTrac WING

Autoclavable KimTrac Wings ensure easy flap retraction with highly improved visibility and accessibility of the operation field.

Item Code: KTW1(Wide type 10EA), KTW2(Long type 10EA), KBJ(KimTrac bending jig)

* KimTrac wings and the KimTrac bending jig are included with the KimTrac.



KimTrac Wing

