




# STS

PROFESSIONAL  
CREATING PRODUCTS.

# HART TEST

HOLLOW AIR REMOVAL TEST FOR SMALL STEAM STERILISERS  
**TO ISO 11140-6:2022**

*STS Professional's  
Hollow Air Removal Test or  
'HART' test (ARTG 460985)  
is manufactured to  
ISO 11140-6:2022 in the  
'helix' style to demonstrate  
the air removal & steam  
penetration efficiency for  
hollow items in small  
steam sterilisers.*



**HART**  
HOLLOW AIR REMOVAL TEST  
**ROUTINE &  
DAILY PCD  
ISO 11140-6**

*The HART test is made up of the  
PCD & 250 Type 2 chemical indicator strips  
& may be used for **daily testing &/or routine  
batch monitoring** in every cycle for hollow  
air removal & steam penetration testing to  
**AS5369:2023.***



REORDER CODE: ST-9S-HART20  
[www.Henryschein.co.nz](http://www.Henryschein.co.nz)