




EASE-IN-SHIELDS™

Protective Loupe Inserts

WORK WITH YOUR LOUPES,
UNINTERRUPTED



You don't need to hold protective glasses in front of your loupes anymore. Ease-In-Shields protective inserts easily slide behind your loupes and fit close to the eyes to provide comfortable, hands-free eye protection. They are available in a variety of options for different laser wavelengths and curing lights.



- ALLOW YOU TO WEAR LOUPES DURING LASER OR CURING PROCEDURES
- PROVIDE HANDS-FREE EYE PROTECTION SO YOU CAN USE BOTH HANDS FOR PROCEDURES
- EASY INSTALLATION: DROP-IN AND GO
- LIGHTWEIGHT POLYCARBONATE CONSTRUCTION



5837 - Dual Wave
OD 5+ @ 800-900nm OD 6+ @ 900-<980nm OD 7+ @ 980-1064nm OD 6+ @ >1064-1080nm



5839 - Curing Light
OD 5+ @ 400 - 500nm



5832 - Universal Soft Tissue
OD 4+ @ 800-<900nm OD 6+ @ 900-<980nm OD 7+ @ 980-1064nm OD 6+ @ >1064-1080nm



5835 - Universal Hard Tissue
OD 4+ @ 1800-3000nm OD 5+ @ 10600nm



5836 - CO₂
0.3<OD<0.49 @ 520-535nm
OD 6+ @ 9000-11000nm



5838 - Diode Soft Tissue
OD 5+ @ 800-830nm



5843 - Pink Diode Soft Tissue Edition
OD 5+ @ 800-830nm



5842 - Multi-Wave
OD 8+ @ 180-385nm OD 1+ @ 650 - <730nm
OD 2+ @ 730 - <800nm OD 4+ @ 800 - <840nm
OD 5+ @ 840 - <960nm OD 7+ @ 960 - 1400nm
OD 4+ @ >1400 - <3000nm OD 6+ @ 10600nm



5844 - Clear CO₂
OD 5+ @ 190 - 390nm
OD 6+ @ 9000 - 11000nm



5845 - Tri-Wave
OD 4+ @ 440-475 nm, OD 2+ @ 655-665 nm,
OD 4+ @ 960-995 nm



5846 - Nd:YAG and Er:YAG
OD 6+ @ 1010-1030 nm, OD 7+ @ 1030-1070 nm,
OD 6+ @ 1070-1080 nm, OD 4+ @ 1800-3000 nm

KIT CONTENTS

- 1 x Ease-In-Shield
- 1 x Adjustable head strap
- 1 x Microfiber cloth
- 1 x Hard case



Call us: AU: 1800.29.09.29 • NZ: 09.801.8123
info.anz@ultradent.com



Ultradent ANZ @ultradentanz

