

KAVO  Kerr

## SonicFill 3



**Product Features:**


**Sonic-activated bulk fill composite that acts as a flowable composite during placement, making the composite easily adaptable to restoration preparations. Rheological modifiers react to the sonic energy applied during placement, reducing viscosity to increase flowability to all cavity surfaces. You will be able to fill most cavities in seconds.**

**The nano-scale filler system (zirconium oxide + silica oxide) will allow you to achieve excellent gloss and wear resistance.**

**Once the sonic activation is removed, the composite reverts back to a viscous, pliable state that is ideal for contouring and polishing.**

**Low shrinkage = void / gap free procedures**

**TRUE one-step bulk fill: 5mm depth of cure**



## Clinician Testimonials





I am interested, what next?

Reach out to your local Kerr Product Specialist Representative to discuss arranging a product demonstration.

Ordering Information

### **36710 SonicFill™ 3 Intro Kit**

The Intro Kit is an ideal way to get introduced to SonicFill 3. It contains a SonicFill Handpiece, a “Getting Started Kit”, a Technique Guide, and a total of 60 Unidose™ tips with 40 Unidose tips in A2, 10 Unidose tips in A1 and 10 Unidose tips in A3.

---

### **SonicFill 3 Unidose Refills**

**Each refill contains 20 Unidose tips each 0.25 g.**

<b>36711</b>	SonicFill 3 Composite A1
<b>36712</b>	SonicFill 3 Composite A2
<b>36713</b>	SonicFill 3 Composite A3
<b>36714</b>	SonicFill 3 Composite B1

---

### **Handpiece and Accessories**

<b>35215</b>	SonicFill Handpiece
<b>35134</b>	Kerr MULTIflex Coupling

[Click for link to product brochure](#)