



Take 1™ Advanced™

VPS Impression Material

Capture outstanding detail in any environment.

Our expertise in high-performance impression materials delivers a VPS material with the optimum combination of physical properties: elasticity, dimensional stability, and the ability to register detail in any environment. Take 1 Advanced has an extensive selection of set times, delivery systems and consistencies on the market.

Features



Outstanding flexibility. Flexible impression material stretches around undercuts without tearing or distorting.



Hydrophilic wash materials with low water contact angle. Displaces intraoral fluids to produce impressions with great detail and no voids or gaps.



Wide selection of viscosities, set times and delivery methods. Helps overcome common difficulties in impression taking, and matches doctor preference.




Excellent dimensional stability. Impression will minimize distortion due to dimensional change or temporary deformation.




Take 1 Advanced Bite Registration

Our Bite Registration provides excellent detail and reduced voids. Take 1 Advanced Bite Registration completes the Take 1 Advanced portfolio providing you a complete solution within a single family of products.



Material Set Times



Data on file.

	LB Wash	Med/Mono	HB Tray	Bite Reg
Options: Delivery, Viscosity & Set Times	A light body material with lowest viscosity and vivid contrasting color.	Medium consistency/ monophase material, a balance between flow and body.	Heavy body material with outstanding dimensional stability and elasticity.	Excellent resistance to deformation to protect occlusal integrity.
Cartridge 50ml 1:1	Fast Set Regular Set	Fast Set Regular Set	Fast Set Regular Set	Work Time: 15 seconds Total Set Time: 1 minute
Volume 380ml 1:1			Fast Set Regular Set	

Ordering Information



Refills

Two Packs: 2 Cartridges (50 ml), 6 Mixing Tips

- 34149 Light Body Wash FS
- 34150 Light Body Wash RS
(Wash viscosity 2-packs also contain intraoral tips)
- 34152 Medium FS
- 34153 Medium RS
- 34158 Tray FS
- 34159 Tray RS
- 34148 Bite Registration

Volume: 2 Hard Cartridges (380 ml)

- 34962 Tray RS
- 34963 Tray FS

Accessories

- 28747 Extruder
- 29363 Mix Tips, Small (pkg of 48)
- 29333 Mix Tips, Large (pkg of 48)
- 27517 Intraoral Syringe Tips, Small Universal (100 disposable)
- 20731 Intraoral Syringe Tips, Large (100 disposable)
- 6000-0000 Sympress Volume Machine
- 34746 Dynamic Mixing Tips for Volume (50 tips, 2 bayonet rings)

Visit kerrdental.com/en-au/ for more information.

Kerr Rep: _____ Email: _____ Phone: _____

Dealer Rep: _____ Email: _____ Phone: _____