





**Diamond | Preparation instruments**  
*Diamante | Instrumentos de preparación*




**Preparation of crown cores**  
 Parallel chamfer  
*Preparación de coronas*  
*Chamfer paralelo*




**KS 1 SC**




		5
<b>Size · Tamaño</b>	∅ 1/10 mm	<b>012</b>
L	mm	7,0

FGSS




● **KS1SC.313. ...** 012

FG





● **KS1SC.314. ...** 012

○<sub>max.</sub> 300000 min<sup>-1</sup>  
 Parallel chamfer  
*Chamfer paralelo*




**KS 2 SC**




		5
<b>Size · Tamaño</b>	∅ 1/10 mm	<b>014</b>
L	mm	8,0

FGSS




● **KS2SC.313. ...** 014

FG





● **KS2SC.314. ...** 014

Parallel chamfer, matches H837KRU carbide prep series  
*Chamfer paralelo, adaptado a la serie de carburo,*  
*H837KRU*




**KS 3 SC**




		5
<b>Size · Tamaño</b>	∅ 1/10 mm	<b>016</b>
L	mm	8,0

FGSS



● **KS3SC.313. ...** 016

FG




● **KS3SC.314. ...** 016


Parallel chamfer, matches H837KRU carbide prep series  
*Chamfer paralelo, adaptado a la serie de carburo,*  
*H837KRU*




**Preparation of crown cores**  
 Occlusal/Lingual reduction  
*Reducción oclusal/lingual*




**KS 5 SC**



		5
<b>Size · Tamaño</b>	∅ 1/10 mm	<b>031</b>
L	mm	6,0

FG



● **KS5SC.314. ...** 031

○<sub>max.</sub> 140000 min<sup>-1</sup>  
 Pear big  
*Pera grande*