


Omni-Matrix[™]

DISPOSABLE RETAINER AND MATRIX











- Unmatched body design is ideal for both single and simultaneous restorations
- Innovative shape allows superior procedural visibility and patient comfort
- Ultra-thin burnishable stainless steel adapts to all preparations
- Unique winged and wingless styles meet individual case needs
- Disposable design saves you time and money

The new Omni-Matrix disposable retainer and matrix is a superior circumferential matrix band solution. It's a simple restorative tool designed to perfectly customize to any preparation. The band's circumference can be easily adjusted simply by twisting the handle and the pivoting head allows it to access any quadrant of the mouth. Once the restoration is complete, the Omni-Matrix band easily releases without disturbing the restorative material.

"Fast, easy, convenient, disposable! What else could you ask for?"
—DR. GEORGE FREEDMAN







WINGED



	Winged 48pk	Wingless 48pk
0.001 Orange	1101	2201
0.0015 Green	1102	2202
Pedo Purple	1104	2204

	Winged 48pk	
Mylar Red	1103	

1. realityesthetics.com.