


Vertex Dental

Baseplates & waxes Vertex™ Light Curing Trayplates



Vertex™ Baseplates



For the production of temporary base plates.

- Thickness 1.3 mm (upper and lower)
- Packaging size: 100 plates per box

Vertex™ Modelling wax



Vertex™ Modelling wax regular is a pink, approximately 1.6 mm thick made of dental wax. These plates have dimensions of about 18.8 x 8.9 cm. For use also in tropical regions we have extra hard wax available as well. Wax Knife also is a part of the Modelling Assortment.

- Bright pink colour
- Packaging: 1000g / 46 wax plates

Vertex™ Preformed wax rims



Preformed wax rims:

- Red hard (100pcs)
- Pink medium (100pcs)

Vertex™ Light Curing Trayplates



Vertex™ Light Curing Trayplates with minimal residual layer after polymerisation makes processing easy. Each box contains 50 trays. Available colours are blue, pink and clear. The Vertex™ Light Curing Trayplates can be polymerised with the Vertex™ Eco Light Box. Also as extra large trayplates available.

The following products are related to the Baseplates & waxes and Trayplates:



Vertex™ Waxknife

Blue waxknife with Vertex inscription.



Vertex™ Eco Light Box

The high intensity and short polymerisation time, makes the machine suitable for polymerisation of individual light curing trays.

Vertex-Dental B.V.

Johan van Oldenbarneveltlaan 62 | 3705 HJ Zeist
P.O. box 10 | 3700 AA Zeist | The Netherlands
Tel. +31 (0)30-69 76 749 | Fax +31 (0)30-69 55 188
info@vertex-dental.com | www.vertex-dental.com